

# ThromboStep

DIAGNOSTIC KITS FOR FLOW CYTOMETRY

CE IVD



Kit designed for quantification of platelet associated immunoglobulin using flow cytometry

## ADVANTAGES

Comparable results with  
**GOLD STANDARD**  
DISCONTINUED KIT

1h25'



High Sensitivity method for Immune Thrombocytopenic screening by flow cytometry.



Differentiate immune from nonimmune thrombocytopenia.



Simple and standardized procedure. Very easy to analyse.



The kit contains all the required reagents to perform the analysis:



2x  
NaHCO<sub>3</sub>



1x  
anti-total  
Rabbit Igs  
FITC



1x  
anti-total  
Igs FITC



1x  
IgA  
FITC



1x  
IgG  
FITC



1x  
IgM  
FITC



1x  
CD42a PE



2x  
Tyrode's  
Buffer

### IMMUNOSTEP SL

Avda. Universidad de Coimbra, s/n  
Cancer Research Center (C.I.C.)  
Campus Miguel de Unamuno  
37007 Salamanca (Spain)  
+34 923 294 827  
info@immunostep.com  
www.immunostep.com

Do you want more information?

Scan this QR code and see  
all the details of our  
ThromboStep Kit.

